MAS® Omni•IMMUNE™ PRO

LIQUID ASSAYED INTEGRATED IMMUNOASSAY CONTROL

IVD

INTENDED USE

Thermo Scientific MAS® Omni•IMMUNE™ PRO is intended for use as an assayed control for monitoring assay conditions in many clinical laboratory determinations. Include **Omni•IMMUNE PRO** with patient serum specimens when assaying for any of the listed constituents. Assay values are provided for the specific systems listed. The user can compare observations with expected ranges as a means of assuring consistent performance of reagent and instrument.

PRODUCT DESCRIPTION

Omni•IMMUNE PRO is a liquid stable control material prepared from human serum. Analyte levels are adjusted with various pure chemicals and preparations from human tissue or body fluids. Prostate specific antigen is from human seminal fluid, alpha-fetoprotein is from human cord serum and carcinoembryonic antigen is from hepatic metastases of human colon adenocarcinoma. Preservatives and stabilizers are added to maintain product integrity.

CAUTION: Omni•IMMUNE PRO is prepared from human source material. Components of the control which are derived from human source material have been tested using FDA accepted methods and found non-reactive for Hepatitis B Surface antigen (HBsAg), Hepatitis C (HCV), HIV-1 and HIV-2. However, no test method can offer complete assurance that products derived from human source material are free of infectious agents. This control must be handled in accordance with recommendations from Centers for Disease Control/National Institutes of Health manual, "Biosafety in Microbiological and Biomedical Laboratories, 2009." The packaging of this product contains dry natural rubber.

DANGER: Omni•IMMUNE PRO control contains ≤0.1% bovine serum albumin (BSA) and ≤0.01% Proclin 300.

H317 - May cause allergic skin reaction.

H334 - May cause allergy or asthma symptoms or breathing difficulties if inhaled.

Avoid breathing mist or vapor. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/eye protection/face protection. In case of inadequate ventilation wear respiratory protection. If on skin: Wash with plenty of soap and water. IF INHALED: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing. If skin irritation or rash occurs: Get medical advice/attention. If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician. Wash contaminated clothing before reuse. Dispose of contents/container to location in accordance with local/regional/national/international regulations.

CONSTITUENT LIST

17-alpha-OH-progesterone	Insulin-like Growth Factor - 1 (IGF-1)
25-Hydroxy Vitamin D	Leuteinizing Hormone (LH)
Acetaminophen	Lidocaine
Adrenocorticotrophic Hormone (ACTH)	Lithium
Alpha-Fetoprotein (AFP)	n-Acetylprocainamide hydrochloride (NAPA)
Aldosterone	Prostatic acid phosphatase (PAP)
Amikacin	PAPP-A*
Anti- Thyroglobulin (A-Tg)	Phenobarbital
Anti-Thyroid Peroxidase (A-TPO)	Phenytion
Benzodiazepine (Nordiazepam)	Phenytion, Free
Beta - Human Chorionic Gonadotropin (β-hCG)	Primidone
Beta-2-Microglobulin	Procainamide
Calcitonin*	Procalcitonin
Cancer Antigen 125 (CA 125)	Progesterone
Cancer Antigen 15-3 (CA 15-3)	Prolactin
Cancer Antigen 19-9 (CA 19-9)	Prostate Specific Antigen (PSA)
Carbamazepine	Prostate Specific Antigen, free (Free PSA)
Carcinoembryonic Antigen (CEA)	Parathyroid Hormone, Intact (PTH, Intact)
Cortisol	Quinidine
C-Peptide	Salicylate
Dehydroepiandrosterone- Sulfate (DHEA-Sulfate)	Sex Hormone Binding Globulin (SHBG)
Digoxin	Triiodothyronine (T3)
Disopyramide	Triiodothyronine, free (Free T3)
Estradiol	Thyroxine (T4)
Estriol, free	Thyroxine, free (Free T4)
Ethosuximide	Thyroxine binding globulin (TBG)
Ferritin	Testosterone
Folate	Theophylline
Fructosamine	Thyroglobulin
Follicle Stimulating Hormone (FSH)	Tobramycin
Gastrin	Tricyclic Antidepressants (TCA) (Nortriptyline)
Gentamicin	Thyroid Stimulating Hormone (TSH)
Growth Hormone	Thyroid Hormone Uptake (T-uptake)
Human Chorionic Gonadotropin (hCG)	Valproic Acid
Homocysteine	Valproic Acid, Free
Immunoglobulin E (IgE)	Vancomycin
Insulin*	Vitamin B12

*Refer to Limitations of Procedures

STORAGE AND STABILITY

Unopened vials of **Omni•IMMUNE PRO** are stable for 30 days from receipt when stored at 2-8°C. Once opened, vials of **Omni•IMMUNE PRO** are stable for 30 days when stored tightly capped at 2-8°C. This product is stable until the expiration date on the box when stored at -25 to -15°C. **Self-defrosting freezers are not suitable.** If any constituent recovery falls outside the published ranges, contact MAS Technical Assistance for current peer group data. Peer comparison data and latest QC lot specific updates are available online through Lab Link xL Quality Assurance Program at www.maslablink.com.

Bacterial contamination produces an increase in turbidity and/or a characteristic odor. Discard vial if evidence of microbial contamination is observed.

CONTROL RANGES

The published control ranges are based upon a combination of replicate assays of representative samples by participating laboratories, instrument/reagent manufacturers and direct correlation with other analytical systems in accordance with established protocol. Instrument values provided are specific to this lot of control only and are intended to assist the laboratory in establishing its own means and ranges. All values have been assigned with instruments and reagents available at the time of assay and expected values may vary with different reagents and/or methodologies. Laboratory established means should fall within the assigned ranges although subsequent instrument, reagent or calibration modifications may invalidate assigned values.

Peer comparison data and latest QC lot specific updates are available online through LabLink® xL Quality Assurance Program at www.maslablink.com. Refer to the Technical Assistance section for contact information.

INSTRUCTIONS FOR USE

Thaw control at room temperature (18-25°C) on a rocker or with periodic gentle inversion until liquid and then immediately store at 2-8°C. Thoroughly mix the contents of the vial before each use by gently inverting for several minutes. Once control is removed from 2-8°C use immediately. Open the vial and transfer the required quantity of control into a clean sample cup. Replace cap immediately and store the opened vial at 2-8°C.

ONCE THAWED, DO NOT REFREEZE THE CONTROL.

Dropper tips are available to help minimize exposure of control samples to air. DO NOT use a syringe needle to withdraw sample through the control bottle cap.

QUALITY CONTROL

All quality control requirements should be performed in conformance with local, state and/or federal regulations or accreditation requirements.

LIMITATIONS OF PROCEDURES

Compatibility of **Omni•IMMUNE PRO** has been demonstrated only with methods shown in this insert. Caution should be employed when using these controls with methods for which values have not been printed.

Insulin, calcitonin and PAPP-A are weighed into **Omni•IMMUNE PRO**; however, no claim is made for expected values nor the stability of these constituents. Testing of and assigning values to the Insulin and PAPP-A constituents of this control is the responsibility of the individual laboratory.

The specificity of antibodies used in immunoassay procedures may vary between lots. As with patient specimens, cross-reactivities between certain constituents in this control serum may occur with some reagents. Consult reagent manufacturer's package insert regarding possible interferences.

Accurate and reproducible results are dependent upon properly functioning instruments, reagents, and good laboratory technique. This product is intended for use as an assayed control for quantitative assays of listed constituents in human serum. This product is not intended for use as a calibrator. For professional use only.

TECHNICAL ASSISTANCE

In the USA, for technical assistance, call 800-232-3342 or 510-979-5417. For insert updates and information, if your laboratory subscribes to LabLink xL, visit www.maslablink.com and select LabLink Extra. Alternatively, to subscribe to LabLink xL call 800-232-3342 or 510-979-5451.

Outside of the USA, if your laboratory subscribes to LabLink xL, visit www.maslablink.com select LabLink Extra. Alternatively, please contact your local sales office or authorized distributor.

Cat. No.	Description	Size
OPRO-101	MAS® Omni•IMMUNE™ PRO, Level 1	6 x 5 mL
OPRO-202	MAS® Omni•IMMUNE™ PRO, Level 2	6 x 5 mL
OPRO-303	MAS® Omni•IMMUNE™ PRO, Level 3	6 x 5 mL
286-606	Dropper Tips	Pkg. 100

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 2024-11-30

thermo
scientific


For Insert updates go to:
www.thermofisher.com/diagnostics

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 2024-11-30

thermo
scientific

For Insert updates go to:
www.thermofisher.com/diagnostics

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
LOT OPRO24111A
OPRO24112A
OPRO24113A

 2024-11-30

For Insert updates go to:
www.thermofisher.com/diagnostics

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
LOT OPRO24111A
OPRO24112A
OPRO24113A

 2024-11-30

For Insert updates go to:
www.thermofisher.com/diagnostics

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
LOT OPRO24111A
OPRO24112A
OPRO24113A

 2024-11-30

For Insert updates go to:
www.thermofisher.com/diagnostics

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
LOT OPRO24111A
OPRO24112A
OPRO24113A

 2024-11-30

For Insert updates go to:
www.thermofisher.com/diagnostics

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

If you have any questions concerning this notification, in the USA please contact the Technical Support Department at 800-232-3342 or 510-979-5417. Outside the USA, please contact your local Subsidiary or Distributor.

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USA
☎ 800-232-3342
sales.diagnostics.fmt@thermofisher.com

Canada
☎ 800-282-4075
info.cddcanada@thermofisher.com

Germany
☎ +49 0800-40 40 771
cdx.de.info@thermofisher.com

Spain, Portugal & Italy
☎ +34 93589 8338
cdx.es.info@thermofisher.com

Asia Pacific
☎ +61 1800 333 110
cdd.asia.info@thermofisher.com

China
☎ +86 800-810-5118
cdx.cn.info@thermofisher.com

New Zealand
☎ 0800 933 966
auinfo@thermofisher.com

Switzerland & Austria
☎ +41 26 663 86 70
cdx.ch.info@thermofisher.com

Australia
☎ +61 1800 333 110
auinfo@thermofisher.com

France
☎ +33 1 40 86 65 20
cdx.fr.info@thermofisher.com

Nordic
☎ +47 2 325 0433
info.nordic.cdd@thermofisher.com

United Kingdom & Ireland
☎ +44 1442 868 940
cdx.uk.info@thermofisher.com

For countries not listed:
distributor.cdd@thermofisher.com

Microgenics Corporation
46500 Kato Road
Fremont, CA 94538 USA
ISO 13485 Certified Company
www.thermofisher.com/diagnostics

B-R-A-H-M-S GmbH
Neuendorfstrasse 25
16761 Hennigsdorf, Germany
Tel: +49 (0) 800 404 077 12
Fax: +49 (0) 800 404 077 13

OPRO-INS-VA
Rev. 10 2020 12

thermo
scientific