




Twelve patients. Twelve viruses.
One testing solution.



xTAG™ Respiratory Viral Panel distributed by Abbott Molecular




THE FIRST AND ONLY FDA-CLEARED VIRAL PANEL
THAT TESTS FOR 12 VIRUSES IN ONE SHIFT

- **Comprehensive** – tests for 12 different viruses and subtypes – covers virtually all major respiratory viruses commonly used in surveillance and clinical settings
- **Sensitive** – excellent sensitivity, including 96.4% for Influenza A, and 100% for RSV A&B
- **Eliminates Guesswork** – aids in ruling in and ruling out viruses suspected in respiratory infection
- **Rapid** – provides results in a single shift
- **Robust** – capable of detecting and identifying multiple pathogens in a sample
- **Cost-effective** – one low price for information on all viral targets

Targets

Viral Family	Viral Subtypes
Respiratory Syncytial Virus (RSV)	A
	B
Influenza A	Non-specific Influenza A
	H1
	H3
Influenza B	
Parainfluenza 1	
Parainfluenza 2	
Parainfluenza 3	
Metapneumovirus	
Rhinovirus	
Adenovirus	

ABBOTT MOLECULAR IS PROUD TO BE A WORLDWIDE DISTRIBUTOR OF xTAG™ RVP. FOCUSED SOLELY ON PROVIDING INNOVATIVE MOLECULAR DIAGNOSTICS THAT SIGNIFICANTLY IMPROVE MEDICAL CARE, ABBOTT MOLECULAR IS A COMPANY THAT YOU CAN TRUST AND HAVE CONFIDENCE IN TO MEET YOUR COMPREHENSIVE RESPIRATORY VIRAL TESTING NEEDS.



For more information about the xTAG™ Respiratory Viral Panel distributed by Abbott Molecular, call 1-800-553-7042 or visit www.abbottmolecular.com

Manufactured for Abbott Molecular by Luminex Molecular Diagnostics
©2008 Luminex Molecular Diagnostics, xTAG, and xMAP are trademarks of Luminex Corporation

